



Rapid Test Range



allied

INTRODUCING

## ALLCARE Rapid Test Range

ALLCARE Rapid test is a chromatographic immunoassay for qualitative detection of specific Antibody and Antigen. Highly Accurate and Easy to use for Quick diagnosis.

### ALLCARE Rapid Test Range



- Dengue NS1+IgG/IgM
- Dengue Ns1
- Dengue IgG/IgM
- Malaria Pf/Pv
- Malaria Pf/Pan
- Troponin I
- HCG Card

### SPECIFICATIONS :

- Quick result within 20 min
- CDSCO Licensed and NABL lab approved
- Control line to ensure validity of assay
- Sensitivity:  $\geq 95\%$  Specificity: 100%
- Shelf life 18-24 Months
- Storage: 4°C – 30°C

**Your partner in quick & reliable diagnosis**

To know more you can reach out us at

C - +91-9699918275

M - info@alliedbiotech.co.in

W - www.alliedbiotech.co.in

BRIDGING THE HEALTHIER FUTURE